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
Table of contents

Introduction	14
Legal Advices	14
Device Description	14
Indications	15
Contraindications	15
Warnings	16
Precautions	16
Application	16
Clinical use	17
Cleaning	17
Sterilization	21
Storage	21
Limited Warranty	22
Disposal	22
Accessories	22

Introduction

We thank you for deciding to purchase AlphaScope™ Hysteroscope of Gimmi® GmbH.

Endoscopes are medical devices that are manufactured to the highest technical standards and require careful handling, care, and storage. Under normal circumstances and with good treatment, they will meet the high demands placed on them and are capable of meeting them for a long time.

Gimmi® GmbH products are precision devices.

Please keep your endoscopes always with the utmost care, so you will enjoy them over a long period.

Read this manual carefully before using your new product. Thereby you keep yourself, the patient and any third parties from damage, which could occur by incorrect installation or by improper operation.

Legal Advices

This technical manual includes ownership protected information, which underlay copyright conditions. All rights are protected. This technical manual must not be copied, neither by photocopies, or other procedures and must not be distributed or saved, neither complete nor in extract and not without explicit written authorization of Gimmi® GmbH. Nomenclature which is at the same time a registered trademark is not especially marked. In the absence of trademark declaration, it may not be gathered that a designation will be a free trademark. Gimmi® GmbH would greatly appreciate being informed about any errors or omissions that may be found in the content of this instruction manual.

Device Description

The AlphaScope™ Hysteroscope is a small diameter, fiberoptic hysteroscope intended for hysteroscopic procedures. The AlphaScope™ Hysteroscope consists of an outer stainless-steel shaft, inner illumination fibers, and an imaging fiber bundle. At the distal end of the imaging bundle is the objective lens and at the proximal end is the eyepiece. Endoscopic light sources are

connected to the light post of the hysteroscope through compatible 3.5 mm light cables. The eyepiece can be connected through a focusing optical coupler to a camera head that carries the image by cable to the camera. The AlphaScope™ Hysteroscope (GMS40A) is intended for use only with the GYNECARE VERSASCOPE™ Sheath (GMS805).

Indications

Hysteroscopes are used as tools for access to the uterine cavity and are not, in and of themselves, a method of surgery.

The AlphaScope™ Hysteroscope is used to permit viewing of the cervical canal and the uterine cavity for the purpose of performing diagnostic and surgical procedures. They are suitable for patients who are deemed suitable by the attending physician, taking into account the indications and contraindications. Generally recognized indications for diagnostic and operative hysteroscopy include:

Diagnostic Hysteroscopy

- Abnormal Uterine Bleeding
- Evaluation of Abnormal Hysterosalpingogram
- Infertility & Pregnancy Wastage
- Amenorrhea
- Pelvic Pain
- Intrauterine Foreign Body

Operative Hysteroscopy

- Directed Biopsy
- Removal of Submucous Fibroids and Large Polyps
- Submucous Myomectomy
- Transection of Intrauterine Septa
- Endometrial Ablation
- Intrauterine Adhesions

Contraindications

Absolute contraindication to hysteroscopy is acute pelvic inflammatory disease (PID). Relative contraindications include:

Hysteroscopy may be contraindicated by the following conditions, depending on their severity or extent:

- Inability to distend the uterus
- Cervical/Vaginal Infection
- Known Pregnancy
- Recent Uterine Perforation
- Cervical Stenosis
- Uterine Bleeding or Menses
- Invasive Carcinoma of the Cervix
- Medical Contraindication or Intolerance to Anesthesia

Contraindications to endometrial ablation

Hysteroscopic endometrial ablation, whether by laser or electrosurgery, should not be undertaken without adequate training, preceptorship, and clinical experience. Additionally, endometrial biopsy should be performed prior to any ablation. The following are clinical conditions that can significantly complicate hysteroscopic endometrial ablation:

- Adenomatous Endometrial Hyperplasia
- Severe Adenomyosis
- Uterine Anomalies
- Uterine Leiomyoma
- Pelvic Pain (Subtle PID)

Contraindications to hysteroscopic myomectomy

Hysteroscopic myomectomy should not be undertaken without adequate training, preceptorship, and clinical experience. The following are clinical conditions that can significantly complicate hysteroscopic myomectomy:

- Severe Anemia
- Inability to circumnavigate a myoma due to myoma size (e.g., predominantly intramural myomas with small submucous components)

Warnings

Federal (USA) law restricts this device to sale by or on the order of a physician.

- Suspicion of pregnancy should suggest a pregnancy test prior to the performance of diagnostic hysteroscopy.
- Only persons having adequate training and familiarity with hysteroscopy should perform hysteroscopic procedures. Consult medical literature relative techniques, complications, and hazards prior to performance of any hysteroscopic procedure.
- When hysteroscopic instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- When using a fluid distension medium, strict fluid intake and output monitoring is required. Excessive intravasation of distension fluid can lead to fluid overload. Potential complications of Continuous Flow Hysteroscopy are:
 - Hyponatremia
 - Hypothermia
 - Pulmonary edema
 - Cerebral edema
 - Uterine perforation resulting in possible injury to bowel, bladder, major blood vessels, and ureter
- Failure to follow all applicable instructions may result in serious surgical consequences.
- Refer to appropriate electrosurgical system user manual indications and instructions to ensure that all safety precautions are taken.
- A thorough understanding of the principles and techniques involved in laser, electrosurgical and ultrasonic procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device and other medical instruments. Ensure that insulation or grounding is not compromised. Do not immerse electrosurgical instruments in liquid unless the instruments are specifically designed and labeled to function in liquid.

Precautions

- Vaginal ultrasonography prior to hysteroscopy may identify clinical conditions that will alter patient management.
- Intrauterine distension can usually be accomplished with pressures in the range of 35–75 mmHg. Unless the systemic blood pressure is excessive, it is seldom necessary to use pressure higher than 75–80 mmHg.

All serious incidents related to the product must be reported immediately to the manufacturer and to the competent authority of the Member State where the user and/or patient is located.

Application


This device is delivered non-sterile and should be cleaned, disinfected, and sterilized immediately prior to the first and each subsequent use according to the procedure provided below.

Using sterile unpacking technique, remove the instrument from the package. To avoid damage, do not flip the instrument into the sterile field. The AlphaScope™ Hysteroscope should be evaluated prior to each use to ensure adequate optical performance. Inspect the image quality using direct visualization or through connection to a video camera. Any signs of foggy or cloudy images indicate that the proximal or distal optics require cleaning.


The instrument should be inspected closely for foreign materials, scratches, cracks, or other signs of damage to the distal tip, rigid sleeve, or eyepiece. If required, the distal lens can be gently cleaned with a cotton swab moistened with isopropyl alcohol.

Only 3.5 mm light cables with Storz® / Olympus® / Gimmi®, Wolf® or ACMI® connectors should be used to connect the hysteroscope to the endoscopic light source.

Use of light cables larger than 3.5 mm could result in extreme heating of the Hysteroscope in the vicinity of the light post. Extreme heating could cause a tissue burn to the user or patient

Do not use the Hysteroscope if there are any signs of damage.

Clinical use

The AlphaScope™ Hysteroscope is introduced into the uterine cavity through the GYNECARE VERSASCOPE Sheath GMS805. Refer to the GYNECARE VERSASCOPE Sheath Instructions for Use for more detailed instructions on insertion.

Cleaning

For reprocessing procedure, the user must take adequate safety precautions (gloves, glasses, etc.).

The AlphaScope™ Hysteroscope must be thoroughly cleaned immediately after each use to remove blood or other foreign material. Dismantle the hysteroscope prior to cleaning. Unscrew the light cable adapters and optical coupler as applicable.

1. General Principle

Note that the manner of treatment can have significant impact on the life of endoscopes. With due care, rigid endoscopes can be used up to 100 times if they are undamaged and can still be cleaned. Any further re-use or the use of damaged and / or contaminated endoscopes is the responsibility of the user.

Check your endoscope optics function (such as adequate lighting of the fibers, clear, sharp, bright, and round picture) and check your endoscope for possible damage (such as sharp edges, loose parts or visible deformation of materials) before each use. However, if you are not sure whether your endoscope optics are operational, we recommend returning the rigid endoscope to our service address for review or replacement. After 100 sterilization cycles, the endoscopes must be sent to our service center for maintenance. Following the maintenance performed by us, the endoscope can be used again for 50 cycles.

2. Instructions for reprocessing

These instructions consider the procedure according to DIN EN ISO 17664 reprocessing of medical devices (cleaning, disinfection, and sterilization). A washer-disinfector (WD) must be used for cleaning and disinfection. It is possible to use a manual cleaning as pre-cleaning process. But because of a significantly lower efficiency and reproducibility of the manual method, it is not allowed to apply a manual pre-cleaning without a subsequent mechanical cleaning and disinfection. Before every cleaning, it is necessary to carry out a pre-treatment. As endoscope optics are provided with irrigation and aspiration stopcocks, these should be disassembled for sterilization. We recommend a steam sterilization / autoclave process. When using a manual cleaning and disinfection procedure, product and process specific validation is the sole responsibility of the user.

3. Preparation before purification

In clinical practice used endoscope optics get sometimes in contact with corrosive etching agents and drugs. Coarse impurities on the instruments must be removed directly after use (within 2 h). Immediately after the operation, a wet cleaning should be performed to prevent the blood, protein, and other substances on the endoscope from drying out and to protect the personnel. Dried protein complicates the cleaning, disinfection, and sterilization. Use running water or a disinfectant solution. The disinfectant should be aldehyde-free (otherwise fixation of blood smears), have proven efficacy (e.g., VAH / DGHM - or FDA approval or CE mark) and should be suitable for disinfecting instruments. To remove impurities manually, use only a soft brush or a clean, soft cloth that is used only for this purpose.

Never use metal brushes, steel wool or other sharp objects!

4. Manual Cleaning/ Disinfection

When choosing the cleaning agent and disinfectant ensure:

- **that these are basically suitable for the cleaning and disinfection of instruments made of metals and plastics,**
- **that a disinfectant with proven efficacy (e.g., VAH/DGHM- or FDA approval or CE mark) is used and that it is compatible with the cleaning agent,**
- **that chemicals used are compatible with the instruments (see "material resistance").**

If possible combined cleaning / disinfecting agents should not be used. Combined cleaning / disinfecting agents can be used only in cases of very low contamination (no visible impurities). Concentrations and contact times of cleaning and disinfectant agents specified by the manufacturer must be strictly adhered to. Only use freshly prepared solutions and sterile or low contaminated water (max. 10 germs/ml), low endotoxin contaminated water (max. 0.25 endotoxin units/ml), e.g. purified water / highly purified water. Additionally, use filtered air for drying.

Cleaning Procedure

- Rinse endoscopes under running water. Remove surface contamination with a soft cloth.
- Rinse the lumens five times with water using a disposable syringe (working channel 50 ml, rinsing and suction channel 10 ml).
- Disassemble the endoscopes as far as possible. Remove the fiber optic connector, adapter, handle, irrigation valves, etc.
- After pre-treatment, place the endoscopes in the cleaning solution for at least 5 minutes so that the instruments are sufficiently covered and thoroughly cleaned with a soft cloth or soft brush. Make sure that the instruments do not touch each other. Move the moving parts back and forth several times during cleaning. Clean the lumens with a soft plastic brush.
- If necessary, connect the rinsing adapter to the endoscopes

- At the beginning and at the end of the cleaning, rinse the lumens five times with cleaning solution (working channel 50 ml, rinsing and suction channel 10 ml).
- Then rinse the products three times with sterile, deionized water.
- Using a disposable syringe, rinse the lumens five times with sterile, deionized water (50 ml working channel, 10 ml irrigation and suction channel).
- For cleaning the endoscopes, use a fresh cleaning solution.
- Check the instruments

Disinfection Procedure

- Disassemble the endoscopes as far as possible. Remove the fiber optic connector, adapter, handle, irrigation valves, etc.
- If necessary, connect the rinsing adapter to the endoscopes.
- Place the endoscopes in the disinfectant solution for 12 minutes so that the instruments are adequately covered. Make sure that the instruments do not touch each other.
- At the beginning and at the end of the disinfection, rinse the lumens five times with disinfectant solution (working channel 50 ml, rinsing and suction channel 10 ml).
- Then rinse endoscopes five times with sterile, deionized water.
- Using a disposable syringe, rinse the lumens five times with sterile, deionized water (50 ml working channel, 10 ml irrigation and suction channel).
- Use a fresh disinfectant solution for the disinfection of the endoscopes
- The evidence of the general suitability of the instruments/devices for effective manual cleaning and disinfection has been provided by an independent accredited test laboratory using the detergent Cidezyme / Enzol and the disinfectant Cidex Opa (Johnson & Johnson GmbH, Norderstedt). The method described above has been considered.

5. Automated cleaning/ Disinfection (WD)

For reprocessing, we recommend our GMS710A screen basket. Pay attention to following points during selection of the WD:

- fundamentally approved efficiency of the WD (for example CE marking according to EN ISO 15883 or DGHM or FDA approval/clearance/registration)
- possibility for an approved program for thermal disinfection (Ao value ≤ 3000 or – in case of older devices - at least 5 min at 90 °C/194 °F; in case of chemical disinfection danger of remnants of the disinfectant on the instruments)
- basic suitability of the program for instruments as well as sufficient rinsing steps in the program
- post-rinsing only with sterile or low contaminated water (max. 10 germs/ml, max. 0.25 endotoxin units/ml), for example purified/highly purified water
- use of filtered air (oil-free, low contamination with microorganisms and particles) for drying
- regularly maintenance and check/calibration of the WD

Pay attention to following points during selection of the cleaning detergent:

- fundamental suitability for the cleaning of instruments made of metallic or plastic material
- additional use - if thermal disinfection is not used - of a suitable disinfectant with approved efficacy (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE mark) compatible to the used cleaning detergent
- compatibility of the used detergents with the instruments (see chapter „material resistance„)

Pay attention to the instructions of the detergent manufacturers regarding concentration, temperature, and soaking time as well as post-rinsing.

Procedure

- Disassemble the endoscopes as far as possible. Remove the fiber optic connector, adapter, handle, irrigation valves, etc.
- If necessary, connect the rinsing adapter to the endoscopes
- Then position the pretreated endoscopes in a washer-disinfector (Disinfector G 7836 CD, Miele, Gütersloh). Make sure that the instruments do not touch each other.
- Connect the lumens with suitable, flexible irrigation tubing and adapters.
- Start the program.
- The machine cleaning was carried out using "neodisher MediClean forte" as a cleaning agent.
- The machine is cleaned at 50 ± 2 °C for at least 5 min. The thermal disinfection was carried out at 90 ± 2 °C for at least 5 min.
- After completion of the machine cleaning / disinfection, remove the endoscopes from the washer-disinfector under low-germ conditions
- The fundamental suitability of the instruments for an effective automated cleaning and disinfection was demonstrated by an independent, governmentally accredited and recognized (§ 15 (5) MPG) test laboratory by application of the WD G 7836 CD, Miele & Cie. GmbH & Co., Gütersloh, (thermal disinfection) and the pre-cleaning and cleaning detergent Neodisher medizym (Dr. Weigert GmbH & Co. KG, Hamburg) considering to the specified procedure.

Failure to promptly clean the hysteroscope after use may result in permanent damage to the instrument and could result in ineffective sterilization.

6. Drying

The endoscope and the camera head must be completely dry after disinfecting. Disinfected endoscopes should always be stored in a closed sterile area or cabinet and protected against heat, radiation, dust, moisture, temperature changes and contamination.

7. Check

Check the device after cleaning or cleaning / disinfection for corrosion, damaged surfaces, impurities, and functionality, respectively. Do not use a damaged device (limiting of the number of re-use cycles see chapter "reusability"). If the device is still dirty, it must be cleaned and disinfected again.

Inspection of the mechanics and endoscope surface

The endoscope surfaces must be undamaged and in particular free of sharp edges. Check for dents, bends, mechanical / thermal damage caused by radiofrequency or laser surgery equipment as well as for cracks and spalling.

Verification of the proximal and distal areas of glass

Glass surfaces must be clean and debris-free. Persistent encrustation, observed during visual examination, should be removed with appropriate cleaning pastes or alcohol-soaked cotton swab or toothpick. Inadequate rinsing of the optics after cleaning and disinfection is often the cause of unusable. Corresponding to the indications appropriate working distance the image must be sharp and clear. A fuzzy, non-circular, cloudy, foggy, image points out to damage.

8. Maintenance

Instrument oils or grease must not be used.

Put the disassembled endoscopes (fibre optic light guide, handles) together again.

9. Packaging

Please insert the cleaned and disinfected device in single-use sterilization packaging's (single or double packaging), which fulfil the following requirements material/process):

- EN ISO/ANSI AAMI ISO 11607 (for USA: FDA clearance)
- suitable for steam sterilization (temperature resistance up to at least 138 °C (280 °F), sufficient steam permeability)
- sufficient protection of the instruments as well as of the sterilization packaging against mechanical damage

Sterilization

Please use only the listed sterilization procedures for sterilization; other sterilization procedures must not be used.

It is user's responsibility to implement the listed sterilization processes to achieve the desired and required sterilization effects.

Steam sterilization:

- fractionated vacuum/dynamic air removal procedure^{1 2} (with sufficient product drying³)
- steam sterilizer according to EN 13060/EN 285 or ANSI AAMI ST79 (for USA: FDA clearance)
- validated according to EN ISO 17665 (valid IQ/OQ (commissioning) and product specific performance qualification (PQ))
- maximum sterilization temperature of 134 °C (273 °F) plus tolerance according to EN ISO 17665
- sterilization time (exposure time at the sterilization temperature):

Area	fractionated vacuum/ dynamic air removal	gravity displacement
USA	at least 4 min at 132 °C (270 °F), drying time at least 20 min	-
other countries	at least 5 min ⁴ at 132 °C (270 °F) / 134 °C (273 °F)	not recommended

It is user's responsibility to validate the process if he wants to use other parameters or substances. Please check with national law.

Storage

Handle with care. This instrument should be stored in a protective container in a dry, clean and dust free environment at modest temperatures of +0°C to 50°C. Protect the distal lens and provided proximal eyepiece from scratches and other physical damage.

¹ at least three vacuum steps

² The less effective gravity displacement procedure must not be used in case of availability of the fractionated vacuum procedure, requires significantly longer sterilization times as well as a sterilizer, procedure, parameter, and product specific validation under sole responsibility of the user.

³ The effectively required drying time depends directly on parameters in sole responsibility of the user (load configuration and density, sterilizer conditions, ...) and is thus to be determined by the user. Nevertheless, drying times less than 20 min must not be applied.

⁴ respectively 18 min (inactivation of prions, not relevant for USA) The flash/immediate use sterilization procedure must not be used. Do not use dry heat sterilization, radiation sterilization, formaldehyde, and ethylene oxide sterilization, as well as plasma sterilization.

Limited Warranty

This device is warranted to the original purchaser against defects in workmanship and material for one year from the date of original purchase to the extent that Gimmi® will, at our option, replace, repair, or refund the purchase price of the device or any part thereof which examination discloses, to our satisfaction, is defective. This warranty does not apply to any product, or part thereof, which has been repaired or altered by persons not authorized by Gimmi®, or which has been subjected to misuse, neglect, or accident. This warranty is In Lieu of all other warranties, express or implied, including, without limitation, any warranties of merchantability, and fitness for a particular purpose and of all other obligations on the part of the manufacturer or seller. We shall not be liable for any anticipated or lost profits, incidental or consequential damages, costs, time charges, or other losses in connection with the device or any of its parts. The use of damaged and/or contaminated Hysteroscope is the responsibility of the user. Disregarding these instructions for use will void the guarantee or warranty claims. We accept no liability in the case of improper handling, incorrect or inadequate preparation or unauthorized repairs. We reserve the right to change or discontinue the equipment at any time without incurring any obligation to make the same or similar changes on equipment previously built or sold by us. Please ensure that instruments returned to Gimmi® for repair or servicing have been thoroughly cleaned, disinfected and sterilized. Gimmi® may refuse to repair uncleaned or contaminated products.

Disposal

Be aware of sharp edges during disposal to prevent injury.

The transposition of European legislation into national law, requires proper disposal of medical devices and instruments. The adjacent symbol means that medical devices must be disposed of separately from household and residual waste. Please dispose of medical products according to local and national laws.

Dispose of contaminated medical equipment at an appropriate collection point. Medical waste is classified according to the Dangerous Goods Regulation for UN number "UN 3291" (Medical Waste)). Please mark the waste container with the adjacent symbol. Dispose of medical waste according to local and national law.

You can dispose of uncontaminated packaging material according to local and national recycling regulations.

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